

For 3 wells: 600 nM, 402 nM, 600 nM, 36 nM, 18 nM, 6 nM
 Protocol for 96-well plate based Oxidative Phosphorylation Inhibitors

Complex I (CI, NADH:Ubiquinone oxidoreductase)

To make from

KCN
 NADH

- 1) Add 2 μ g mitochondria to x wells of a 96 well plate and dilute to 140 μ L. Include at least three positive control wells and three negative control wells which include rotenone.
- 2) To three of the six wells add 20 μ L of potassium phosphate buffer (0.5 M, pH 7.5), and 12 μ L of fatty acid-free BSA (50 mg/ml), and 6 μ L of KCN (10 mM) and 2 μ L of NADH (10 mM) total mic add 40
- 3) Adjust the volume to 198.8 μ L with distilled water. - 198.8 - 180 = 18.8 m
- 4) In three control wells add in parallel the same reagents from steps 2-3 but with 2 μ L, 1 mM rotenone. Add rotenone wells first and wait 30s
- 5) In parallel add the same reagents from steps 2-3 but include 2 μ L of test compounds at 10⁻³ three-fold dilutions. Double the amount of NADH then add normal wells
 - a. Concentrations should be calculated using 200 μ L final volume. I.e. to get 10 μ M final concentration, make a 1 mM stock solution.
- 6) Mix plate on shaker and read baseline at 340 nm for 2 min. Specific complex I activity is the rotenone-sensitive activity.

Complex II (CII, succinate dehydrogenase)

DCPIP

use 200
 Mitos / reaction

- 1) Add 120 μ L ddH₂O, 10 μ L potassium phosphate buffer (0.5 M, pH 7.5), 4 μ L of fatty acid-free BSA (50 mg/ml), 6 μ L KCN (10 mM), 5 μ L of succinate (400 mM), and 0.5 μ g of mitochondria, and 29 μ L of DCPIP. for 20 wells - 2400 ddH₂O, 200 pot, 80 m BSA, 120 KCN, 100 Succ
- 2) Adjust the volume to 198.8 μ L with ddH₂O. 10 m mto
- 3) Mix plate in instrument using shaker for 8 minutes and read for 2 min at 600 nm.
- 4) Start the reaction by adding 0.8 μ L of 12.5 mM DUB, mix and follow the decrease of absorbance at 600 nm for 3 min.
- 5) Check the specificity of complex II activity by running the assay after the addition of 10 μ L of either 1 M malonate or 50 mMTTFA before starting the reaction. The degree of inhibition usually exceeds 85% with TTFA and 95% with malonate, indicating a high specificity of reaction.

Complex III (decylubiquinol cytochrome c oxidoreductase)

KCN

use 1200
 well

- 1) Add 146 μ L of ddH₂O, 10 μ L of potassium phosphate buffer (0.5 M, pH 7.5), 15 μ L of oxidized cytochrome c, 10 μ L of KCN (10 mM), 4 μ L of EDTA (5 mM, pH 7.5), 2 μ L of Tween-20 (2.5% vol/vol) and 0.1 μ g mitochondria to x wells of a 96 well plate and dilute to 140 μ L. 182 total
- 2) Prepare, in parallel, at least three control wells containing the same reagents as above, but with the addition of 2 μ L of a 1 mg/ml antimycin A stock.
- 3) Adjust the volume to 198 μ L using ddH₂O and mix by shaking.
- 4) Read baseline at 550 nm for 2 min.
- 5) Start the reaction by adding 2 μ L of a 10 mM decylubiquinol solution, mix rapidly by and immediately read the increase in absorbance at 550 nm for 1-2 minutes. Specific complex II activity is the antimycin A-sensitive activity.

Complex IV (CIV, cytochrome c oxidase)

KCN

- 1) Add 80 μ L of ddH₂O, 100 μ L of phosphate buffer (0.1 M, pH 7.0), 12 μ L of reduced cytochrome C (1 mM) and read the baseline activity at 550 nm for 2 minutes.

for 35 wells

42.8 m ddH₂O, 5 m pot, 760 m reduced cy c

for 15 wells 1200, 1500, 180

- 2) Adjust the volume to 199 μL
 - a. 7 μL ddH₂O
- 3) Start the reaction by adding 1 μL of sample (0.1 μg of mitochondria), mix and monitor the decrease of absorbance at 550 nm. Should take 3 minutes.
- 4) To check the specificity of complex IV, add 30 μL of 10 mM KCN in a separate reaction

4.1.5 μM

For 3.487 $\mu\text{M}/\text{mL}$

↳ need 15 μM

$$\text{↳ } 0.1 \text{ nM} \times 15 = 1.5 \text{ nM} / 3.487 \text{ nM/mL}$$

= 0.44 μL into

14.56 μL total